

croquis

SYB-MD-CRO-202511V1

croquis®

Samyang Biopharm Corporation

| | |
|----------------------|---|
| Business Areas | 4 |
| Major Achievements | |
| Major Location | 5 |

croquis®

| | |
|----------------------|----|
| Brand Motive | 8 |
| MONO | 10 |
| SCREW | 12 |
| VOLUME | 14 |
| BARBED CUTTING | 16 |
| BARBED MOLDING | 18 |

Samyang Aesthetics


| | |
|----------------------|---|
| Art Technology | 6 |
|----------------------|---|

Guaranteed Quality

| | |
|-------------------------------|----|
| Storage Test | 20 |
| BHS Test | 21 |
| * BHS (Barb Holding Strength) | |

Brand Vision

| | |
|---|----|
| SMART Academy | 22 |
| (Medical Aesthetics Education Platform) | |




Samyang Biopharm Corporation

Through accumulated experience and technology

With more than 25 years of experience in suture manufacturing and polymer technology, Samyang Biopharm now intends to expand its portfolio for Medical Aesthetics.

Business Areas



Major Achievements



Major Location



Samyang Discovery Center
(Headquarter & R&D Center)

| | |
|--------------------|--|
| Completion | 2016 |
| Related Department | Marketing & Sales / Regulatory Affairs / Planning / R&D Center |
| Location | Seongnam-si, Gyeonggi-do, Korea |

Pharmaceuticals Plant



COMPLETION
1997

PRODUCT
Injection, API

LOCATION
Daejeon, Korea

CERTIFICATIONS
EU GMP, PMDA (Japan), GMP, KGMP (Korea)

Medical Devices Plant



COMPLETION
1996

PRODUCT
Sutures, Hemostats, Patches

LOCATION
Daejeon, Korea

CERTIFICATIONS
KGMP, CE Mark, ISO13485 (Medical Devices), KGMP (Pharmaceuticals)

Samyang Biopharm USA Inc.




COMPLETION
2018

BUSINESS MODEL

- NRDO (No Research Development Only)
- Therapeutics for Oncology and Rare Diseases
- License-in / Out

LOCATION
Cambridge, MA, USA




Samyang Aesthetics

Everyone is born beautiful.

The definition of beauty has changed from wanting to look like someone to bringing out your beauty within. Samyang believes in value as art. Samyang Aesthetics is derived from combining art with technology. By bringing out your embedded beauty with Samyang's technology, Samyang Aesthetics provides a helpful and safe solution to stay young and beautiful.


Art Technology




croquis®

Brand Motive

croquis® is motivated by the drawing technique of describing an object in a short period of time. Find your embedded beauty through skin rejuvenation and tissue repositioning.



RETURNING TO THAT YOUTH



Tightening ●●● | Volume ●○○○

croquis MONO stimulates neocollagenesis, which not only improves your overall skin quality but also fully absorbs into the body over time.

Material : PDO (Polydioxanone)
Contains : 60ea / box (10ea / Al pouch * 6)
Thread Type : _____

| Needle | | Thread | |
|--------|-------------|--------|-------------|
| Gauge | Length (mm) | USP | Length (mm) |
| 25G | 60 | 5 - 0 | 90 |
| | 90 | | 150 |
| 27G | 60 | 6 - 0 | 90 |
| | 40 | | 50 |
| 29G | 50 | | 70 |
| | 25 | | 30 |
| 30G | 25 | | 30 |
| | 40 | | 50 |





Tightening ●●● | Volume ●●○

croquis[®] SCREW stimulates neocollagenesis and improves your overall skin and wrinkle through volumization effect. SCREW type also gets absorbed completely into the body over time.

Material : PDO (Polydioxanone)
 Contains : 60ea / box (10ea / Al pouch * 6)
 Thread Type : 

| Needle | | Thread | |
|--------|-------------|--------|-------------|
| Gauge | Length (mm) | USP | Length (mm) |
| 26G | 50 | 5 - 0 | 70 |
| | 60 | | 90 |
| | 90 | | 130 |
| 27G | 40 | 6 - 0 | 52 |
| | 50 | | 70 |
| 29G | 40 | 6 - 0 | 52 |
| | 50 | | 70 |

B RINGING THAT SMOOTH AND TIGHT SKIN





croquis[®] VOLUME effectively adds volume to deep wrinkles or fold, providing a smoother and more rejuvenated appearance.



Material : PDO (Polydioxanone)

Contains : 30ea / box (6ea / AI pouch * 5)



Thread Type :

| L-Cannula | | Thread | |
|-----------|-------------|-------------|-------------|
| Gauge | Length (mm) | USP | Length (mm) |
| 19G | 40 | 6 - 0 x 9ea | 60 |
| | 60 | | 90 |
| 21G | 40 | 7 - 0 x 9ea | 60 |
| | 60 | | 90 |

A DDING THAT VOLUME TO Y OUR SKIN



PULLING AND HOLDING THAT BEAUTY



BARBED CUTTING

Skin Quality ●●● | Lifting ●●●

croquis[®] BARBED CUTTING is used for tissue fixation and repositioning. With absorbable component and neocollagenesis, BARBED CUTTING type provides a younger and beautiful look.

Material : PDO (Polydioxanone)
Contains : 30ea / box (6ea / Al pouch * 5)

Thread Type :  
(3D, Bi-Directional) (Reverse Bi-Directional)

| L-Cannula | | Thread | | Barb |
|-----------|-------------|--------|-------------|--------------------------|
| Gauge | Length (mm) | USP | Length (mm) | |
| 18G | 100 | 2 | 150 | (3D, Bi-Directional) |
| 19G | 60 | 0 | 105 | |
| | 90 | | 140 | |
| 20G | 60 | 2 - 0 | 105 | (Reverse Bi-Directional) |
| | 90 | | 140 | |
| 19G | 40 | 2 | 50 | |
| | 50 | | 70 | |

 **BARBED MOLDING**

Skin Quality ●●● | Lifting ●●●+

croquis[®] BARBED MOLDING is used for a stronger lift. By using the pine tree-like barb shape, BARBED MOLDING type gives a helpful solution to skin laxity.



Material : PDO (Polydioxanone)

Contains : 30ea / box (6ea / Al pouch * 5)

Thread Type : 

| L-Cannula | | Thread | | Barb |
|-----------|-------------|--------|-------------|---|
| Gauge | Length (mm) | USP | Length (mm) | |
| 18G | 100 | 2 | 150 |  (2D, Bi-Directional) |
| 19G | 100 | 1 | 185 | |
| 19G | 70 | 1 | 155 |  FIX (2D, Multi-Directional) |
| | 100 | | 185 | |
| 20G | 70 | 0 | 155 | |
| | 100 | | 185 | |

EN HANCING AND KEEPING THAT LIFT



Guaranteed Quality

High stability of croquis[®] can solve quality-related issues during storage and distribution processes. Such stability provides mutual trust to global customers.

1) Storage Test

Result after 5-year accelerated aging test with competing brands


MONO
SCREW
TYPE : MONO
USP : 6 - 0



BARBED CUTTING
USP : 3 - 0, 4 - 0




BARBED MOLDING
USP : 2 - 0



2) BHS Test

croquis[®] BARBED CUTTING type's high BHS demonstrates effective facial lifting.


BARBED CUTTING
BHS : Barb Holding Strength (kgf)
USP : 0



Brand Vision

croquis® as Market Standard

Samyang Aesthetics provides high and safe quality products through well-controlled and strict manufacturing process.




Education Platform

SMART Academy

SMART Academy is an online training platform for medical aesthetics experts.

- Product Information
- Guide Protocol
- e-Course: Procedural Video



www.SMART-ACADEMY.co.kr

Samyang Medical Aesthetics Research & Training

Draw your beauty,
Capture the moment

croquis®



croquis[®]

samyang



Samyang Aesthetics



SMART Academy



E-brochure

